



# Elvaplain

PLAIN CATGUT

Natural | Absorbable | Monofilament

KEMENKES RI AKD 21603323578

# Elva Surgical Suture

## Elvachromic CHROMIC CATGUT

## Elvaplain PLAIN CATGUT

## Elvasilk SILK

Specification	Monofilament, Natural	Monofilament, Natural	Multifilament, Natural
Color	Brown	Beige	Black
Complete Absorption	70-90 days	56-63 days	Nonabsorbable
Tensile Strength (Post Implantation)	2 weeks, approx 50% 3 weeks, approx 30%	1 weeks, approx 50%	Indefinite
USP Size	USP 1 - USP 4-0	USP 1 - USP 4-0	USP 1 - USP 4-0

Benang Bedah **Elvaplain** adalah benang bedah natural yang dapat diserap, monofilamen, yang terbuat dari jaringan ikat yang dimurnikan (kolagen) yang berasal dari lapisan serosa sapi (*bovine*) atau lapisan fibrosa submukosa usus domba (*ovine*). Benang Bedah **Elvaplain** dikemas dengan larutan yang mengandung isopropanol dan air.

Benang Bedah **Elvaplain** diindikasikan untuk digunakan dalam aproksimasi jaringan lunak secara umum dan/atau ligasi, termasuk penggunaan dalam prosedur mata tetapi tidak untuk prosedur kardiovaskular dan neurologis.

Benang bedah **Elvaplain** yang didistribusikan, berwarna dalam ukuran USP 1 sampai USP 4-0 dalam panjang yang bervariasi, dengan jarum.

Produk ini dipasarkan dan didistribusikan di bawah entitas **PT Enseval Medika Prima**.

**Elvaplain** Surgical Suture is an absorbable, monofilament, natural surgical suture made of purified connective tissue (collagen) derived from either the serosal layer of beef (*bovine*) or the submucosal fibrous layer of sheep (*ovine*) intestines. **Elvaplain** Surgical Suture is packaged with a solution containing isopropanol and water.

**Elvaplain** Surgical Suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures but not for use in cardiovascular and neurological procedures.

**Elvaplain** surgical sutures are supplied in USP sized 1 through 4-0 in a variety of lengths, with needle.

This product is marketed and distributed under the entity **PT Enseval Medika Prima**.



# DESKRIPSI UMUM KODE BENANG BEDAH

## GENERAL DESCRIPTION OF SURGICAL SUTURE CODE

**Contoh** Example





### 1 Inisial Material *Material Initial*

H = Chromic Catgut

### 2 Ukuran Benang *Thread size*

Kode	Ukuran
9	USP 2
1	USP 1
0	USP 0
2	USP 2-0
3	USP 3-0
4	USP 4-0

### 3 Lengkungan *Curvature*


Kode	Lengkungan
1	1/2 
2	3/8 

### 4 Panjang Jarum *Needle Length*










Kode = Panjang Jarum  
Code = Needle Length








### 5 Bentuk Ujung Jarum *Point Shape*






### 6 Nomor Berjalan *Running Number*

- 10** : Dus isi 12, Panjang Benang 75 cm  
Box of contents 12, Suture Length 75 cm
- 11** : Dus isi 18, Panjang Benang 75 cm  
Box of contents 12, Suture Length 75 cm
- 12** : Dus isi 24, Panjang Benang 75 cm  
Box of contents 12, Suture Length 75 cm

Principal Code	EMP Code	Box of Content	USP Size	Needle Length	Needle Description
L1148P-00	L148P	12	1	 48 mm	1/2 Taper Point
L1148P-01		18			
L1148P-02		24			
L1140P-10	L140P	12	1	 40 mm	1/2 Taper Point
L1140P-11		18			
L1140P-12		24			
L0140P-10	L040P	12	0	 40 mm	1/2 Taper Point
L0140P-11		18			
L0140P-12		24			
L0136P-00	L036P	12	0	 36 mm	1/2 Taper Point
L0136P-01		18			
L0136P-02		24			
L2125P-10	L225P	12	2-0	 25 mm	1/2 Taper Point
L2125P-11		18			
L2125P-12		24			
L2126P-10	L226P	12	2-0	 26 mm	1/2 Taper Point
L2126P-11		18			
L2126P-12		24			
L2135T-00	L235T	12	2-0	 22 mm	1/2 Taper Cutting
L2135T-01		18			
L2135T-02		24			

Principal Code	EMP Code	Box of Content	USP Size	Needle Length	Needle Description
L2137P-10	L237P	12	2-0	 37 mm	1/2 Taper Point
L2137P-11		18			
L2137P-12		24			
L2140P-10	L240P	12	2-0	 40 mm	1/2 Taper Point
L2140P-11		18			
L2140P-12		24			
L2224R-10	L224R	12	2-0	 24 mm	3/8 Reverse Cutting
L2224R-11		18			
L2224R-12		24			
L3125P-10	L325P	12	3-0	 25 mm	1/2 Taper Point
L3125P-11		18			
L3125P-12		24			
L3225R-10	L325R	12	3-0	 25 mm	3/8 Reverse Cutting
L3225R-11		18			
L3225R-12		24			
L3119P-10	L319P	12	3-0	 19 mm	1/2 Taper Point
L3119P-11		18			
L3119P-12		24			
L3122T-10	L322T	12	3-0	 22 mm	1/2 Taper Cutting
L3122T-11		18			
L3122T-12		24			

Principal Code	EMP Code	Box of Content	USP Size	Needle Length	Needle Description
L4117P-10	L417P	12	4-0	 17 mm	1/2 Taper Point
L4117P-11		18			
L4117P-12		24			
L4216R-10	L416R	12	4-0	 16 mm	3/8 Reverse Cutting
L4216R-11		18			
L4216R-12		24			
L4219R-10	L419R	12	4-0	 19 mm	3/8 Reverse Cutting
L4219R-11		18			
L4219R-12		24			



Manufactured by :

**PT Forsta Kalmedic Global**

A Kalbe Company

Jl. Pulogadung No.23, Industrial Estate (JIEP), Jakarta 13930

Phone : 021 - 5086 7788

<https://forsta.co.id>



Distributed by :

**PT Enseval Medika Prima**

A Kalbe Company

Jl. Pulo Lentut No.12, Industrial Estate (JIEP), Jakarta 13930

Phone : 021 - 4682 3234

<https://emp.co.id>